OPERATING INSTRUCTIONS

EN

1. Twenty to thirty minutes before use, remove the product box from the refrigerator, remove the blister pack from the box and allow the syringe to come to room temperature. Do not use if the sterile package appears to be damaged.

2. **Peel off the blister backing.** (Fig. 1)

3. While holding the blister open side down, bend the blister and allow the syringe to fall gently onto the sterile surface.

Alternatively, hold the blister open side up and bend back the blister until the barrel's luer end is exposed. Gripping the luer end of the barrel, remove the syringe from the blister.
Do not remove the syringe from the plunger end. (Fig. 2)

4. **To remove the tip cap from the syringe:**

- Hold the syringe upright on the ribbed part of the white closure.

- With the other hand, take hold of the white cap of the closure system. Tilt back and forth carefully until the cap disconnects and can be pulled off after the seal is broken. (Fig. 3.1)

Do not rotate.

- Remove the cap in a straight upward direction. Do not touch the syringe tip. (Fig. 3.2)

- Do not use the syringe if the closure cap perforation is already broken.

5. **Attach the cannula.** (Fig. 4.1, Fig. 4.2)

Attention: Do not apply pressure to the plunger rod while the cannula is being affixed. Verify that the cannula is properly locked to the Luer Lock Adaptor (LLA). Do not overtighten the LLA; this can lead to loosening of the LLA from the barrel.

6. **Apply gentle pressure to the plunger** in order to expel air from the cannula and verify that the syringe is operating properly. (Fig. 5)

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Each ml of BioLon® Prime contains:	
Sodium hyaluronate	12 mg
Sodium chloride	8.5 mg
Disodium hydrogen phosphate dodecahydrate	0.56 mg
Sodium dihydrogen phosphate dihydrate	0.05 mg
Water for injections	q.s

DESCRIPTION
BioLon Prime is a sterile 1.2% solution of highly purified, high molecular weight viscoelastic sodium hyaluronate in phosphate-buffered saline.
Sodium hyaluronate is a physiological material found in both animal and human tissues.
BioLon Prime is a very highly purified product extracted from bacterial cells. It is a polysaccharide consisting of a repeating disaccharide of N-acetylglucosamine and sodium glucuronate, linked by alternating β-1,3 and β-1,4 glycosidic bonds.
Chemically identical in all species, hyaluronate can be found in the vitreous and aqueous humor of the eye, in synovial fluid, in skin and in the umbilical cord.

When introduced into the anterior segment of the eye during surgical procedures, BioLon Prime maintains a deep anterior chamber throughout surgery, facilitates surgery and reduces trauma to the corneal endothelium and surrounding tissues. Its viscoelasticity also helps to repel the vitreous face and discourage formation of a postoperative flat chamber.

Any traces of BioLon Prime left in the anterior segment of the eye after surgery dissipate through Schlemm’s canal within a week.

INDICATION
Surgery in the anterior segment of the eye.

CONTRAINDICATIONS
When used as recommended, there are no known contraindications.

PRECAUTIONS
Precautions are limited to those normally associated with the surgical procedures being performed.

Instilling excessive amounts of BioLon Prime into the anterior or posterior segments of the eye may increase intraocular pressure. Increased postoperative intraocular pressure may also be caused by a pre existing glaucoma condition or by compromised outflow and by operative procedures and sequelae thereto, including enzymatic zonulysis, absence of an iridectomy, trauma to filtration structures, and by blood and lenticular remnants in the anterior chamber.
Because these factors vary from case to case and are difficult to predict, the following precautions are recommended:

Do not overfill the eye with BioLon Prime.

At the end of surgery, all remaining BioLon Prime should be removed by irrigation or aspiration.

Carefully monitor the intraocular pressure, especially during the immediate postoperative period.
If a significant increase in pressure is observed, treat appropriately.

Avoid the trapping of air bubbles behind BioLon Prime.

Since the presence of minute quantities of impurities (proteins, etc.) cannot be totally excluded, the physician should be aware of immunological, allergic and other potential risks associated with the injection of biological substances.

On rare occasions, viscoelastic products containing sodium hyaluronate have been observed to become slightly opaque or to form a slight precipitate upon instillation into the eye. The clinical significance, if any, of this phenomenon is not known.
The physician should, however, be aware of this possibility, and, should it be observed, the cloudy or precipitated material should be removed by irrigation and/or aspiration.
Avoid re-use of the cannula.
If re-use is necessary, rinse the cannula thoroughly to remove all traces of residual materials.

DOSAGE
Dosage varies with type of surgery.
Usually a dose of 0.2 to 0.7 ml is injected into the anterior segment of the eye.
Greater amounts are used in the posterior segment.

INCOMPATIBILITIES
Mixing of quaternary ammon